

LATEX FREE | SINGLE PATIENT USE | DISPOSABLE

HORIZON™

HYPERINFLATION

Horizon™ Hyperinflation System

Deliver a high standard of care while effectively managing healthcare costs with the single use Horizon Hyperinflation System. Each system comes standard with an integrated pressure manometer. Optional components and comfort features allow you create a custom resuscitator to meet the needs of each patient and of your facility's budget.


Hyperinflation System

- Systems equipped with a stay-put dial to enable clinician to achieve targeted static pressure
- Three pleated breathing bag sizes available: 1/2, 1 and 2 liter
- Non-inflatable face masks available



Hyperinflation System

Series	Description	PK
1000HZ	Hyperinflation System Standard	10
2000HZ	Hyperinflation System + Manometer + Pop-Off	10
3000HZ	Hyperinflation System + Manometer	10



Integrated Pressure Manometer
helps promote proper pressure delivery

2000HZ Series
shown with integrated pressure manometer and pop-off valve

Patient Port Swivel
allows clinician to easily move around patient while maintaining proper resuscitation



Stay-Put Dial
enables clinician to achieve targeted static pressure

Convenience Components



Face Masks
Non-inflatable, air-filled, soft cushion masks are available in variety of sizes

Safety Components



Integrated Pressure Manometer
The innovative color-coded dial clearly denotes pressure delivery:
Green Target pressure level <20cmH₂O
Yellow Potentially dangerous pressure level
Red Dangerous pressure level



Pop-Off Valves
40cm/H₂O & 25cm/H₂O pop-off provides an extra safety feature to prevent peak inspiratory pressures from reaching excessive levels further protecting the infants fragile lungs

 **SOVEREIGN MEDICAL INC.**

919-644-1113 orders@sovmed.com www.sovmed.com