

Distributed by



SOVEREIGN MEDICAL, INC

919-644-1113

sovmed.com



eitan
medical
sapphire

sapphire™ FastTest PM™

Your Solution for Independent
Onsite Infusion Pump
Maintenance.



30 Tests
per kit





Want to streamline your facility's maintenance process for infusion pumps and avoid relying on outside service centers?

FasTest PM™ is the solution

FasTest PM™ is an on-site maintenance servicing tool for the Sapphire Infusion Pump Platform, which enables you to:



Quickly and easily conduct full maintenance verification tests, including occlusion detection, air detection, and pump accuracy



Reduce maintenance costs and equipment downtime, which puts you in control



Receive an immediate certification report helping maintain compliance

FasTest PM gives you independence and control of the maintenance process

Standardize pump testing with a dedicated tool and software. FasTest PM offers you the option to do annual PM tests between treating patients, keeping more pumps available and verified for patient use.

FasTest PM makes pump maintenance easier than ever

Utilize this easy and intuitive solution, which provides step-by-step, guided navigation, allowing your staff to easily and quickly learn the software.

FasTest PM brings efficiency home

Save costs, time, and improve efficiency by moving your annual PM tests in-house. No more chasing down pumps. Within minutes, you can have your manufacturer certification, ensuring your infusion devices are verified.

FasTest PM helps you streamline your maintenance processes

Deploy an efficient maintenance system that is fast, cost-effective and sustainable. FasTest is reusable and offers multiple applications, giving your team full independence to conduct maintenance on your Sapphire fleet.



Technical features

Sequence of tests

Visual inspection, Safety features testing, Calibration verification

Average PM run time

< 5 minutes

Maximum simultaneous PM processes

8

Calibration verification process

- Air Detection
- Flow Accuracy
- Occlusion Sensing

Kit content

- Administration set - 1800 mm tubing, with: spike, Syringe, back check valve, Y-connector, T-connector, administration cassette, slide clamp, Sleeve, and Fluid bag cap
- 500 mL fluid bag

FasTest PM certificate

- Automatically generated certificate
- Device SW version
- Certification kit S/N
- Test date
- Technician name
- Visual inspection & safety features test results
- Air Detection, accuracy & occlusion test results



SOVEREIGN MEDICAL, INC

620H Valley Forge Rd - Hillsborough, NC 27278

919-644-1113 sovmed.com

Technical Service Contact Information:

service@eitanmedical.com

Customer Support Contact information:

support@eitanmedical.com

+972.73.2388888 | fax: +972.73.2388800

www.eitanmedical.com