

NeoTree™ Tube Support



The **NeoTree™ Tube Support** works in conjunction with neonatal corrugated breathing tubes. Reversible to secure tubes for term and pre-term neonates and supports tubes in all positions. Assists in the reduction of tension on the endotracheal tube or face

mask by helping to prevent kinks in the tube, accidental extubation and/or movement of the endotracheal tube. Made of lightweight aluminum with a baked-on powder coat finish, durable and can be easily cleaned.

Product Code: 420-1948

Distributed by

SOVEREIGN MEDICAL, INC

620H Valley Forge Rd - Hillsborough, NC 27278

919-644-1113

sovmed.com



**Reversible to
Secure Term and
Pre-term Neonatal
Corrugated
Breathing Tubes**

**Lightweight
Aluminum**

**Easy to Clean
Powder Coat
Finish**

Reusable

Durable