

FLO₂MAX Filtered Oxygen Mask

from Solutions in Critical Care

Closed Mask Breathing System with 99.99% Filtration Efficiency




The FLO₂MAX Filtered Oxygen Mask is an all-in-one closed mask breathing system that offers low, medium and high oxygen concentration and nebulization from just one mask. All patient exhalation is filtered through a submicron, hydrophobic 3M™ filter that maintains respiratory isolation during oxygen and drug aerosol therapy to help protect healthcare workers and patients from airborne infections. The FLO₂MAX is a cost-effective, high-quality solution to oxygen therapy.

Highlights

- Deliver 30% - 99% Oxygen
- Compatible with standard nebulizer for drug therapy
- Over 99.99% efficiency (BFE/VFE)
- Three mask sizes to ensure best fit possible
- 100% Latex-Free

How It Works



Item ID	Description
6100	Large Adult FLO ₂ MAX Mask with O ₂ Reservoir Tubing
6105	Small Adult FLO ₂ MAX Mask with O ₂ Reservoir and Tubing
6110	Child FLO ₂ MAX Mask with O ₂ Reservoir and Tubing
6250-N	Rescuer Nebulizer



Solutions in Critical Care

Distributed by

SOVEREIGN MEDICAL, INC

620H Valley Forge Rd

Hillsborough, NC 27278

919-644-1113 sovmed.com

Now on GPO Contracts

Premier Contract Number: PP-DS-743

Vizient Contract Number: DM0170