


FLUSSO™

TFI

To Enable a Safer Circuit Disconnect



Made in Canada



MMSI

McARTHUR
MEDICAL SALES INC.



Flusso™ TFI allows for temporary interruption of gas flow to maintain a patient's positive end expiratory pressure (PEEP) and minimize potential exhaled aerosols during ventilator circuit disconnect.



Flusso™ Highlights:

- Lightweight
- Clear housing permits visualization of function
- Low deadspace (1.86 cc)
- Compact
- Minimal effort to activate
- Safety tabs to prevent over compression

Product Code 420-1979 25/case

FLUSSO™
TFI

Manufactured By:

McArthur Medical Sales Inc
1846 5th Concession Rd W
Rockton, ON LOR 1X0
www.mcarthurmedical.com

MMSI
McARTHUR
MEDICAL SALES INC.

Distributed by

SOVEREIGN MEDICAL, INC

620H Valley Forge Rd
Hillsborough, NC 27278

919-644-1113 sovmed.com

Flusso™ is a trademark used under licence by McArthur Medical Sales Inc. US and International Patents Pending

TFI-030521Rev1