

Pedi-Fit™

NON-INVASIVE PATIENT INTERFACE

Right-Fit Pediatric Patients with Excellent Seal and Comfort

The Artemis Pedi-Fit Non-Invasive Pediatric Nasal and Full-Face Masks were specifically designed to right-fit respiratory patients in the Pediatric Intensive Care Unit (PICU). The soft, molded silicone material and flexible forehead pad ensure excellent seal and patient comfort during CPAP and VPAP therapy.

Pedi-Fit Features

- Molded silicone mask cushion provides excellent seal and comfort
- Unique clip pivots the mask to conform to the patient’s nose, preventing air from leaking into the patient’s eyes
- Flexible, one-piece forehead pad evenly conforms to the forehead for stability
- Quiet air vent minimizes noises and enables efficient CO₂ washout
- Each design is available in four sizes enabling healthcare providers to right-fit each patient

Item ID	Description	U/M
ART00PF	Pedi-Fit Nasal Mask, Small (21 mm x 18.6 mm)	Each
ART0PF	Pedi-Fit Nasal Mask, Medium (22 mm x 24.10 mm)	Each
ART01PF	Pedi-Fit Nasal Mask, Large (26.5 mm x 28.2 mm)	Each
ART02PF	Pedi-Fit Nasal Mask, Adult X-Small (30.5 mm x 32.3 mm)	Each
ART200PF	Pedi-Fit Full-Face Mask, Small (28mm x 37 mm)	Each
ART20PF	Pedi-Fit Full-Face Mask, Medium (35 mm x 43 mm)	Each
ART21PF	Pedi-Fit Full-Face Mask, Large (41 mm x 50 mm)	Each
ART22PF	Pedi-Fit Full-Face Mask, X-Large (46 mm x 60 mm)	Each
TST600	72” Ventilator Circuit with Smooth Bore Tubing	20/Bx



Nasal Mask



Full-Face Mask



Ventilator Circuit

Distributed by

SOVEREIGN MEDICAL, INC

620H Valley Forge Rd - Hillsborough, NC 27278
919-644-1113 www.sovmed.com



Solutions in Critical Care
(866) FOR-SCC1 or (866) 367-7221
www.solutionsincriticalcare.com


NOW ON GPO CONTRACTS

Premier Contract Number: PP-DS-743
Vizient Contract Number: DM0170


Pedi-Fit™

NON-INVASIVE PATIENT INTERFACE

Sizing Chart for Pedi-Fit Nasal Mask



Sizing Chart for Pedi-Fit Full-Face Mask



Distributed by

SOVEREIGN MEDICAL, INC

620H Valley Forge Rd - Hillsborough, NC 27278
919-644-1113 www.sovmed.com



Solutions in Critical Care

(866) FOR-SCC1 or (866) 367-7221
www.solutionsincriticalcare.com

NOW ON GPO CONTRACTS

Premier Contract Number: PP-DS-743

Vizient Contract Number: DM0170