

**Solutions in  
Critical Care**

# CPAP


**with Inspiratory/Expiratory HEPA Filter**



## Provide Maximum Patient Benefits & Infection Control with One Device

The Solutions in Critical Care (SCC) CPAP System was designed to deliver optimal  $FiO_2$  levels while conserving oxygen resources. Known as the disposable CPAP with the lowest oxygen consumption, the SCC CPAP is also the only disposable CPAP with an inspiratory/expiratory filter for infection control. With multi-function access ports, the SCC CPAP is cost-efficient, easy-to-use, delivers maximum patient results and minimizes the risk of infection for patients and healthcare providers.

### SCC CPAP System Performance Data



## SCC CPAP Features

- Lowest Oxygen Consumption of Any Disposable CPAP: 5-10 Lpm
- Filter captures 99.997% of 1  $\mu$ m particles (TB and Bacteria) and 99.993% of .28  $\mu$ m particles (Viruses) that are inhaled and exhaled.
- Port for Oxygen Tubing Connection to Boost  $FiO_2$
- Port for Nebulizer Tubing Connection to Deliver Medications
- Calibrated Pressure Manometer
- Low Cost to Buy and Operate

### SCC CPAP Inspiratory/Expiratory Filter

#### Performance Data


**99.997%**

1  $\mu$ m particles (TB/Bacteria)

**99.993%**

.28  $\mu$ m particles (Viruses)

Distributed by

 **SOVEREIGN MEDICAL, INC**  
620H Valley Forge Rd - Hillsborough, NC 27278  
919-644-1113 sovmed.com



**Solutions in Critical Care**

40 Boroline Road, Allendale, New Jersey 07401  
www.solutionsincriticalcare.com • (201) 512-1100

**NOW ON GPO CONTRACTS**

Premier Contract Number: PP-DS-743  
Vizient Contract Number: DM0170