

from Solutions in Critical Care

# Ventilator Circuit

## 72" Disposable Ventilator Circuit with Smooth Bore Tubing

A high-quality, low-cost ventilation circuit compatible with adult and pediatric Bi-Level and CPAP systems.



### Technical Specifications

- 1 each 180 cm smooth bore tubing - ID:19.5mm(30.5mm) OD:23.5mm(30.5mm)
- 2 each PVC straight connector - 19M-22F
- 1 each Straight Connector with tethered Sampling Port - 22M/15F-22M/15M
- 1 each 2m sampling line - ID:3.5mm OD:5.3mm; Wall Thickness:1.8mm
- 1 each hose hanger (certain configurations only)

Item ID	Description	U/M
TST600	72" Ventilator Circuit with Smooth Bore Tubing	20/Bx

## Contact Your Local Distributor

Distributed by

### SOVEREIGN MEDICAL, INC

620H Valley Forge Rd - Hillsborough, NC 27278

919-644-1113 [www.sovmed.com](http://www.sovmed.com)



**Solutions in Critical Care**

(866) FOR-SCC1 or (866) 367-7221  
[www.solutionsincriticalcare.com](http://www.solutionsincriticalcare.com)

**NOW ON GPO CONTRACTS**

Premier Contract Number: PP-DS-743

Vizient Contract Number: DM0170