

eitan
medical
sapphire™

Sapphire™

Infusion Solutions Catalog



Distributed by:

SOVEREIGN MEDICAL, INC

620H Valley Forge Rd - Hillsborough, NC 27278

919-644-1113 sovmed.com



Eitan Medical p. 3

Infusion kits p. 4

Administration sets p. 5

Accessories p. 8

Spare and maintenance parts p. 10

Service training and maintenance tools p. 11

Clinical training p. 12

Care is the starting point of our innovation at Eitan Medical, at the heart of the development of our drug administration and infusion solutions.

We develop innovative technologies for the benefit of our patients, as well as our partners, healthcare providers, caregivers and technicians.

Eitan Medical is committed to supporting patients around the world and improving care by promoting comfort and independence throughout the treatment cycle.

With a global network specialising in drug delivery and infusion solutions, Eitan Medical is here to support you, our customers, our partners and the patients we serve together.

Sapphire, one infusion pump, multiple therapies.



Hospitals



Home care



Infusion centers



Emergency services



Military



Infusion pump kits



Sapphire™ Multi-Therapy infusion pump kit

Included in the kit



Sapphire™ Multi-Therapy Pump



Mini cradle with knob key



Power supply

P/N 17000-028-0072



Sapphire™ Epidural infusion pump kit

Included in the kit



Sapphire™ Epidural pump



Bolus handle



Mini cradle with splitter and knob key



Power supply

P/N 17000-031-0035

Administration sets



General use

NAME	PRODUCT NUMBER	DESCRIPTION	LENGTH (in)	PRIME VOLUME	UNITS PER BOX
AP403-01	12003-000-0013	Microbore set with a vented/non-vented spike, Fx administration cassette, slide clamp, and spin male luer lock.	114	~6 mL	30
AP409-01	12005-000-0003	Macrobore set with a vented/non-vented drip chamber with 15µ filter, back check valve, needleless Y-site, roller clamp, Fx administration cassette, slide clamp, needleless Y-site, and spin male luer lock.	109	~23 mL	30
AP433-01	12000-000-0044	Macrobore set with a vented/non-vented drip chamber with 15µ filter, Fx administration cassette, roller clamp, needleless Y-site, and spin male luer lock.	113	~22 mL	20
AP423-01	12003-000-0023	Microbore set with a non-vented spike, Fx administration cassette, slide clamp, back check valve, and spin male luer lock.	115	~6 mL	30

Filter sets

NAME	PRODUCT NUMBER	DESCRIPTION	LENGTH (in)	PRIME VOLUME	UNITS PER BOX
AP223-01	12003-000-0033	Microbore administration set with vented/non-vented spike, Fx administration cassette, back check valve, 0.2µ filter, slide clamp, needleless y-site and spin male luer lock.	117	~9 mL	30
AP204-01	12003-000-0070	Microbore set with a non-vented spike, Fx administration cassette, 0.2µ filter, slide clamp, back check valve, and spin male luer lock.	116	~8 mL	30

Administration sets



Filter sets

NAME	PRODUCT NUMBER	DESCRIPTION	LENGTH (in)	PRIME VOLUME	UNITS PER BOX
AP240-01	12003-000-0069	Microbore set with a non-vented spike, Fx administration cassette, 1.2µ filter, slide clamp, back check valve, and spin male luer lock.	116	~8 mL	20
AP 206 -01	12000-000-0022	Macrobore set with a non-vented spike, Fx administration cassette, 1.2µ filter, slide clamp, and spin male luer lock.	116	~23 mL	20

Specialty Sets

NAME	PRODUCT NUMBER	DESCRIPTION	LENGTH (in)	PRIME VOLUME	UNITS PER BOX
AP424-01	12003-000-0024	Microbore administration set with vented/non-vented spike, Fx administration cassette, slide clamp, back check valve, y-site and spin male luer lock, and with female luer lock, slide clamp and a back check valve.	116	~7 mL	30
AP402-01	12004-000-0004	Microbore yellow marked set with vented/non-vented spike, Fx administration cassette, slide clamp and spin male luer lock.	114	~6 mL	30
AP214-01	12004-000-0002	Microbore yellow marked set with a vented/non vented spike, Fx administration cassette, 0.2µ filter, slide clamp, and spin male luer lock.	116	~8 mL	30

Administration sets





Specialty Sets

NAME	PRODUCT NUMBER	DESCRIPTION	LENGTH (in)	PRIME VOLUME	UNITS PER BOX
AP417-01	12004-000-0022	Microbore yellow marked set with a vented/non vented spike, Fx administration cassette, and neuraxial (NRFit®) connector.	114	~6 mL	30
AP133-01	12006-000-0010	Macrobore blood set with two non-vented spikes, two roller clamps, drip chamber with 200 µ filter, Fx administration cassette, roller clamp, needleless Y-site, and spin male luer lock.	115	~30 mL	20
AP504-01	12000-000-0045	Macrobore UV resistant set with an opaque vented/non-vented drip chamber with 15 µ filter, roller clamp, Fx administration cassette, needleless Y-site, and spin male luer lock.	113	~26 mL	20
AP416-01	12003-000-0012	Microbore set with a female luer lock, Fx administration cassette, slide clamp, needleless Y-site, and spin male luer lock.	93	~5 mL	30
AP410-01	12005-000-0007	Macrobore set with a vented/non-vented drip chamber with 15 µ filter, roller clamp, spin male luer lock, and with an infusion bag hanger.	28.1	~9 mL	30
AP411-01	12003-000-0078	Macrobore administration set with vented/non-vented spike, slide clamp, needleless y-site, Fx administration cassette, slide clamp, needleless y-site and spin male luer lock.	72	~13 mL	30



Support

<p>Mini cradle with splitter</p>  <p>P/N 17000-026-0007</p>	<p>Mini cradle without splitter</p>  <p>P/N 15042-000-0010</p>	<p>Multi pump mounting system</p>  <p>P/N 17000-020-0008</p>	<p>Splitter</p>  <p>P/N 15043-000-0005</p>
---	--	---	--

Power supply

<p>Power supply</p>  <p>P/N 15072-000-0007</p>	<p>External battery pack</p>  <p>P/N 15079-111-0001</p>	<p>Integrated power supply (IPS)*</p>  <p>P/N 15072-000-0010</p>
<p>Sapphire power cord compatible with IPS*</p>  <p>P/N 15070-112-0005</p>	<p>Sapphire Multi-pump power cord compatible with IPS</p>  <p>P/N 15171-112-0002</p>	

* Order both PNs for a complete IPS solution







Bags

<p>Infusion pouch 500 mL</p>  <p>P/N 05080-605-0001</p>	<p>Travel case</p>  <p>P/N 15020-010-0002</p>	<p>Homecare large backpack</p>  <p>P/N 15088-000-0001</p>
---	---	---

Pain management



<p>PCA lockbox 250 mL</p>  <p>P/N 15063-000-0001</p>	<p>PCA lockbox 250 mL strap and pouch</p>  <p>P/N 15063-000-0004</p>	<p>PCA lockbox 100 mL</p>  <p>P/N 15061-131-0004</p>	<p>Bolus handle</p>  <p>P/N 15071-000-0005</p>
--	--	--	--



<p>Li-Ion battery</p>  <p>P/N 05020-160-0001</p>	<p>Connection knob flange Key for mini cradle</p>  <p>P/N 05040-380-0019</p>	<p>Connection knob key for mini cradle</p>  <p>P/N 05040-380-0018</p>	<p>Connection knob tip bolt for mini cradle</p>  <p>P/N 05040-380-0027</p>
--	--	--	--



<p>Press plane for mini cradle</p>  <p>P/N 05040-380-0026</p>	<p>100 mL lockbox key</p>  <p>P/N 05070-190-0015</p>	<p>250 mL lockbox key</p>  <p>P/N 10010-122-0046</p>
---	--	---

Maintenance Accessories

<p>Sapphire FasTest PM Kit (for 30 tests)</p>  <p>P/N 16024-000-0002</p>	<p>Communication cable</p>  <p>P/N 05020-110-0213</p>
--	---



Additional products

<p>Carrying case</p>  <p>This is not an Eitan Medical product</p> <p>P/N 200513-1480N5-A</p>	<p>RS232 to USB 4 port adapter</p>  <p>P/N 15077-102-0001</p>
--	---

Pump management & maintenance tools

The Sapphire future-ready, software-based platform expands beyond the pump to provide a complete solution for pump management and maintenance with a variety of software-based, automated time-saving tools.

 <p>FasTest PM*</p> <p>Perform comprehensive maintenance verification tests quickly and easily.</p>	 <p>Sapphire Configuration Manager (SCM)*</p> <p>Configure a pump and duplicate its parameters across your entire fleet.</p>	 <p>Sapphire Pump Loader *</p> <p>Access the latest developments in Sapphire software.</p>
 <p>Event Viewer *</p> <p>View and download up to 8,000 events stored in log memory.</p>	 <p>Sapphire Drug Library</p> <p>Create and edit a drug library of up to 12,000 drugs, in up to 40 different Clinical Care Areas</p>	 <p>Bill of Health **</p> <p>Diagnose and test the basic functionality of your pumps.</p>
 <p>Clean Up Tool *</p> <p>Reset the software of your devices to a basic setting.</p> <p>*Level 1, **Level 2</p>		



Eitan Medical addresses your educational needs

A customized learning experience based on
your educational needs



Virtual and
in-person
instructor-led
training



Interactive learning
experience with Eitan
Medical Academy



Variety of resources
with video tutorials
and user guides

Train your team, **Track** your progress, **Provide** safe patient care

Lined area for taking notes, consisting of multiple horizontal dotted lines.

10050-003-0009-MMTL Rev.08/01.2023

Distributed by

SOVEREIGN MEDICAL, INC

620H Valley Forge Rd - Hillsborough, NC 28278
919-644-1113 sovmed.com



www.eitanmedical.com